

Corsinel[®]
by Tytex



HIGHLY
BREATHABLE



4-WAY
STRETCH



SEAMLESS
KNITTING



LATEX
FREE

STOMASAFE PLUS

OSTOMY SUPPORT BELT



ABOUT STOMASAFE PLUS

- ✕ **DOUBLE LAYERED DEVICE** WITH A HOLE FOR THE OSTOMY POUCH
- ✕ **OUTER LAYER RETAINS THE OSTOMY POUCH** SAFELY, PREVENTING THE POUCH FROM MOVING
- ✕ **HELPS MUFFLE OR REDUCE NOISE** FROM STOMA BAG

OEKO
TEX[®]

STANDARD
100

6234-7001
DTI

EUROPEAN-BASED PRODUCTION
FOR OPTIMAL QUALITY,
SUPPLY CHAIN EFFICIENCY
AND FLEXIBILITY



Corsinel®

by Tytex

STOMASAFE PLUS
ARTICLE NUMBER:3216




AVAILABLE COLOURS



COMPOSITION

91% PA
9% EA

CARE HANDLING



- Wash with similar colours and in a washing bag
- Use of chlorine based bleach or chlorine based disinfection will decrease durability and function



AVAILABLE IN SIZES



ABOUT THE PRODUCTS

The StomaSafe Plus Ostomy Support Belt conceals and supports the ostomy pouch, enabling the user to live an active life. It helps to reduce content noise and can be used with both 1- and 2-piece systems.

FEATURES

StomaSafe Plus is not just a belt; it is a dependable daily companion, offering vital support to secure the ostomy pouch. With StomaSafe Plus, the wearer can confidently navigate their day with a profound sense of security.

Developed with a keen appreciation for discretion, StomaSafe Plus excels at both fixating and concealing the ostomy pouch. Its high-quality construction ensures a soft, skin-friendly experience, boasting a seamless design for optimal comfort.

Ease of use is paramount with StomaSafe Plus. Featuring a soft, stretchy fabric, it is effortless to put on, ensuring a hassle-free experience every time. This flexible material is not only highly breathable but also exceptionally comfortable, free from irritants.

Prioritizing health, StomaSafe Plus is latex-free, safeguarding the well-being of the wearer. It also minimizes noise from the pouch and its contents, providing a discreet and peaceful experience. Moreover, StomaSafe Plus seamlessly accommodates both 1- and 2-piece ostomy systems, demonstrating its versatility and adaptability to unique requirements.

The commitment to safety is unwavering. StomaSafe Plus has undergone rigorous testing to meet the stringent standards set by Oeko-Tex standard 100, assuring the highest quality and peace of mind. StomaSafe Plus is the natural choice when wanting a secure, comfortable, and discreet ostomy support solution.

INDICATION

For post-surgical use after construction of a colostomy, an ileostomy or a urostomy.
For use with both 1- and 2-piece systems

PRECAUTION

- Follow the Instructions for Use
- Follow recommendations given by your doctor or ostomy nurse
- Use the device according to the indications to prevent injury or damage
- The device may influence the stoma drainage if it is too small or fitted too tight
- Remove the device if any unexpected discomfort, stoma or skin irritation occurs
- By urgent need for removal of the device, open with scissors
- No chemicals added during treatment/finishing of the device.
- The device may contain softening and antistatic agents up to 1.5 %

RECOMMENDED SIZES

SIZE RANGE	S/M		M/L		L/XL	
	CM	IN	CM	IN	CM	IN
StomaSafe Plus	85 - 110	33 - 43	105 - 125	41 - 49	120 - 140	47 - 55



Measurement at the widest abdominal circumference

WWW.CORSINEL.COM

CE TYTEX A/S HAS BEEN CERTIFIED ACCORDING TO ISO 9001, ISO 14001 AND ISO 13485