



WE CAN

Provide Options for More
Comfort and Flexibility



Soft Convex Wafer Depth = 4mm (1/8")

ESTEEM™ Soft Convex

Flexibility combined with shallow convexity
— designed to deliver comfort.

Features include:

- **Flexible Convexity:** conforms to the body and helps improve the pouch seal by minimizing leaks
- **Easy-View Window:** for accurate positioning and observation
- **Belt Tabs:** attach an ostomy support belt for added security
- **Filter:** may help minimize odors and ballooning*
- **Safe Seal Clipless Closure:** may help make emptying and cleaning easy*
- **Anti-Reflux Valve:** may help prevent backflow**

*Drainable options only **Urostomy options only



Convex



Designed to be
Skin Friendly



Flexible,
Conforms to Body



Designed to
Minimize Pressure

ESTEEM™ Soft Convex

Esteem™+ Soft Convex combines flexibility with a shallow convexity, and is designed to conform well to the body and be comfortable to wear.



The Esteem™+ Soft Convex system is available in 3 plateau widths to address a range of stoma management issues.



V-Size Plateau Guidance

Convexity Pressure

A narrow plateau width may be used to create pressure directly around the stoma.

Plateau Width

A wider plateau width may be used to help flatten out skin creases away from the stoma.



Esteem™+ Soft Convex | Drainable

(Pouch length = 8 2/3 in)

HCPCS: A4389



Drainable | Cut-to-Fit

Product Code	Plateau Size	Cutting Range	Quantity
422363	V1	20-47mm (3/4 - 1 7/8 in)	10
422364	V2	15-40mm (3/5 - 1 3/5 in)	10
422365	V3	10-28mm (3/8 - 1 1/8 in)	10

Esteem™+ Soft Convex | Urostomy

(Pouch length = 7 2/3 in)

HCPCS: A4392



Urostomy | Cut-to-Fit

Product Code	Plateau Size	Cutting Range	Quantity
422551	V1	20-47mm (3/4 - 1 7/8 in)	10
422557	V2	15-40mm (3/5 - 1 3/5 in)	10
422562	V3	10-28mm (3/8 - 1 1/8 in)	10



Drainable | Pre-Cut

Product Code	Plateau Size	Hole Diameter	Quantity
422539	V1	25mm (1 in)	10
422540	V1	30mm (1 1/8 in)	10
422541	V1	35mm (1 3/8 in)	10
422545	V2	25mm (1 in)	10
422546	V2	30mm (1 1/8 in)	10



Urostomy | Pre-Cut

Product Code	Plateau Size	Hole Diameter	Quantity
422552	V1	25mm (1 in)	10
422553	V1	30mm (1 1/8 in)	10
422558	V2	20mm (3/5 in)	10
422559	V2	25mm (1 in)	10
422560	V2	30mm (1 1/8 in)	10
422564	V3	20mm (3/5 in)	10
422565	V3	25mm (1 in)	10



Enrolling your patients is simple. Contact your me+™ concierge at 1-877-585-0470 or visit mepluscare.com